# 7.15 Determination of serum/ plasma urea

## Principle

A pink colour complex is formed when urea is treated with di-acetyl monoxime and thiosemicarbazide in the presence of sulphuric acid, phosphoric acid and ferric chloride which is read at 525 nm in spectrophotometer.

### Reagents

- Reagent A: 300 ml concentrated H<sub>2</sub>SO<sub>4</sub> (95-98%) 50 ml distilled water 100 ml concentrated H<sub>3</sub>PO<sub>4</sub> (85%) 100 mg ferric chloride (FeCl<sub>3</sub>) Mix well and volume it to 1 litre with distilled water.
- Reagent B: 500 g diacetyl monoxime (DAMO) 10 mg thiosemicarbazide (TSC) Mix well with distilled water to make 100 ml solution



Fig. 7.13 Centrifuging of extracts

- 3. Reagent C: Prepare immediately before use with Reagents A & B at 2:1 proportion.
- 4. Trichloro acetic acid (TCA) 5% (w/v) with distilled water.
- 5. Standard urea solution (0-150 nmol) Working standard (50 μg/ml): Dissolve 50 mg urea in 1 litre of distilled water.

#### Procedure

- 1. Protein is precipitated out from serum/plasma 0.2 ml mixed with 1.8 ml of 5 per cent TCA.
- 2. Centrifuge at 2000 rpm for 10 min (Fig. 7.13) and collect supernatant for determination using the protocol given below.
- 3. Boil in water bath for 5 minutes.
- 4. Cool to room temperature and read absorbance at 525 nm against reagent blank.
- 5. Plot a standard curve and read concentration of the sample against absorbance.

#### Protocol

Reagents/	Blank	Sample	Standard				
Solutions (ml)			1	2	3	4	5
Protein free sample	0.20	-	-	-	-	-	-
Standard	-	-	0.02	0.04	0.06	0.08	0.10
Distilled water	-	0.20	0.18	0.16	0.14	0.12	0.10
Reagent C	3.00	3.00	3.00	3.00	3.00	3.00	3.00

#### Calculation

OD unknown

Urea (mg/ 100 ml) = ----- x concentration of standard x 50 OD standard

**Reference:** Rahmatullah, M. and Boyde, T.R.C. 1980. An improvement in determination of urea using diacetylmonoxime method with and without deproteinization. *Cli. Chem. Acta,*. 107: 3-9.